

2022



ANESTHESIA

Endotracheal Tube (Oral/Nasal)



- Ref no:1.3.1
- Is inserted into the trachea for the primary purpose of establishing and maintaining a patent airway and to ensure the adequate exchange of oxygen and carbon dioxide
- Available both with cuff and without cuff.
- Made of transparent thermosensitive PVC that softens at body temperature.
- With high-volume, low-pressure cuff.
- Packing:1pc/blister bag,100pcs/ctn

	Endotrachealtubus
	Tube endotrachéal
	Tubo endotraqueal
	気管内チューブ

Reinforced Endotracheal Tube (Oral/Nasal)



- Ref no:1.3.2
- Available both with cuff and without cuff.
- Made of transparent thermosensitive PVC that softens at body temperature.
- With high-volume, low-pressure cuff.
- A stainless steel spring is inserted in the tube to minimize the risk of kinking or crushing.
- Packing:1pc/blister bag,100pcs/ctn

	Verstärkte Endotrachealtubus
	Tube endotrachéal renforcé
	Reinforced Endotracheal Tube
	補強気管内チューブ

Tracheostomy Tube



- Ref no:1.3.3
- Printed neck plate with size information for easy identification.
- Straps provided in the pack for fixation of the tube.
- Smooth rounded tip of the Obturator reduces trauma during insertion.
- High volume, low pressure Cuff provides excellent sealing.
- Rigid blister pack provides maximum protection for the tube.
- 1pc/bag,50pcs/ctn or as your required.

	Verstärkte Endotrachealtubus
	Tube endotrachéal renforcé
	Reinforced Endotracheal Tube
	補強気管内チューブ

Disposable PVC Laryngeal Mask



- Ref no:1.3.4
- The tube is made of medical grade PVC.
- The cuff is made of medical grade silicone.
- Suitable for adult, child and infant use.
- 1 pc/bag, 50 pcs/ctn or as required.

	Larynxmaske
	Masque laryngé
	Máscara laríngea
	喉頭マスク

Disposable Silicon Laryngeal Mask



- Ref no:1.3.5
- Made of medical grade Silicon.
- Suitable for adult, child and infant use.
- 1 pc/bag, 50 pcs/ctn or as your required.

	Larynxmaske
	Masque laryngé
	Máscara laríngea
	喉頭マスク

Reusable Silicone Laryngeal Mask



- Ref no:1.3.6
- Made of medical grade Silicon. Suitable for adult, child and infant use.
- Autoclave sterilization only, can be reused up to 40 times.
- 1 pc/bag, 50 pcs/ctn or as your required.

	Larynxmaske
	Masque laryngé
	Máscara laríngea
	喉頭マスク

Endobronchial Tube



- Ref no:1.3.7
- Left side,right side
- Consists of 2 suction tubes,2 angular connectors,1 Y connector, Double lumen endobronchial tube
- Packing:1sets/blister,30sets/ctn

	Endobronchiale Rohr
	Fin b Ronchia l Tsube
	Ende b Ronchia l Tsube
	気管支内チューブ

Sengstaken-Blakemore Stomach Tube



- Ref no:1.3.9
- Nature latex material,
- 3Way
- Size:18
- Sterile
- Packing:1pc/blister,20pcs/.ctn

	Magensonde
	Estomac Tube
	Tubo de estómago
	胃管

Endotracheal Tube Holder



- Ref no:1.3.8
- For Comfortable and Secure Stabilization
- This product should be used by a specialist. After the endotracheal tube is placed into trachea , the tube is fixed with this holder. Firstly put endotracheal tube into the center hole of theholder, retain the tube in the appropriate position, then use the screw to fasten the tube. Finally fix the holder around a patient's chin and the base of the patient's head with strap.
- Single use only.
- Packing:1pc/blister,100pcs/ctn

	r an der t Rachea l Tsube ho l Ende
	Tube Holder endotrachéale
	r al final t Rachea l Tsube ho l
	気管内チューブホルダー

Anesthesia Mask



- Ref no:1.3.10
- Size:0#-5#
- Administer anesthetic gases to a patient through inhalation
- With favourable transparency.
- It is 100% latex free, complying with biocompatibility standards.
- Soft air cushion ensures flexibility, air tightness and low odor.
- Packing: 1pc/polybag, 50pcs/ctn or as your required

	Anästhesie-Maske
	Masque d'anesthésie
	Masque d'anesthésie
	麻酔マスク

Oral Pharyngeal Airway Berman Type



- Ref no:1.3.12
- Size:40-120mm
- A traumatic soft-rounded edges
- Two sides channel design permits passage of suction catheters to clear secretions from the back of the throat without obstructing the air passage.
- Combined with 8 different sizes in a transparent PE bags
- Material:High density polyethylene
- Packing:1pc/polybag, 500pcs/ctn

	Berman Atemwegs
	Voies respiratoires berman
	Berman vías respiratorias
	バーマン気道

Guedel Type Airway




- Ref no:1.3.11
- Size:40mm-110mm
- Used to maintain a patent (open) airway. It does this by preventing the tongue from covering the epiglottis, which could prevent the person from breathing. When a person becomes unconscious, the muscles in their jaw relax and allow the tongue to obstruct the airway-Color coded for easy identification of sizes.
- Made of polypropylene
- Packing:1pc/polybag, 500pcs/ctn

	Guedel
	Guedel voies respiratoires
	Cánula de Guedel
	Guedel タイプ気道

B.V. Filter




- Ref. no.1.3.13.1
- Breathing filter with heat&&moisture exchanger are used at the patient end connection of a breathing system to prevent heat and moisture loss. The filters are designed to prevent bacterial/viral cross contamination via breathing system used on an infectious patient during ventilation.
- 75*70mm, PP material, Filtration membrane
- HV Filtration efficiency: 99.999%
- Packing: 1pc/bag, 200pcs/ctn

	B.V. Filter
	B.V. filtre
	Filtro B.V.
	フィルタB.V.

B.V. Filter



- Ref. no.1.3.13.2
- Bacterial&viral's retention; protect the breath between patients and the machine.
- Resistance (Pa): 50 at 30 l/min
- Bacterial Retention: 99.99%
- Viral Retention: 99.999%
- Packing: 1pc/bag, 200pcs/ctn

	B.V. Filter
	B.V. filtre
	Filtro B.V.
	フィルタB.V.

Closed Suction Catheter



- Ref no:1.3.14
- Size: FR6, FR8, FR10, FR12, FR14, FR16, FR18
- 24hours
- Designed for removing bronchial secretions in intubated Patients during mechanical ventilation.
- Designed for 24 hours use
- Better control of infections

	Absaugkatheter
	Cathéter d'aspiration
	Catéter de succión
	吸引カテーテル

Disposable Anesthesia Breathing Circuit



- Ref no:1.3.15
- A closed system interface delivers oxygen and anesthetic gas to patient.
- Economical dual limb circuit for single patient use.
- Lightweight for patient's comfort.
- Collapsible corrugated tubing for more flexibility and easy storage.
- Length can be made as your request
- Packing: 1pc/polybag,50pcs/ctn or as your required

	Atemschlauch
	Circuit respiratoire
	Circuitos
	呼吸回路

CPAP Nasal Mask



- Ref no:1.3.16
- Can provide continuous positive airway pressure.
- Can provide bi-level positive airway pressure treatment for adult patients.
- For sleep apnea.
- For respiratory failure caused by various reasons.
- For chronic obstructive pulmonary disease nocturnal hypoxemia

	CPAP-Maske
	CPAP MASK
	La máscara de CPAP
	CPAP のマスク

CPAP Full Face Mask







- Ref no:1.3.17
- Size:S,M,L
- The ultimate in economy and comfort for CPAP therapy
- Conforms comfortably to the user's face without the need for forehead pad
- Allows the user the ability to read, watch TV, and feel less confined

	CPAP-Maske
	CPAP MASK
	La máscara de CPAP
	CPAP のマスク

Catheter Mount Tube with Angle Piece



- Ref no:1.3.18
- clear plastic, corrugated, compressible, light weight tube, 90 mm length, with 22mm female termination and angle piece with 22F/15F termination
- Packing:1pc/blister,200pcs/ctn

	Catheter Berg
	Catheter Mont
	Catéter Monte
	カテーテルマウント

Intubating Stylet



- Ref no:1.3.19
- Made of malleable aluminum covered with plastic sheath.
- Soft rounded tip for minimizing patient trauma during intubation
- Flexible use for assistance of intubation
- Packing:1pc/blister bag,100pcs/ctn

	Intubating Stilett
	Stilet intubation
	Estilete de intubación
	挿管スタイルット

Endotracheal Tube Introducer (Bouge)



- Ref no:1.3.20
- Smooth elbow design of tip can assist intubation, decreasing the damage of soft tissue.
- Scale mark indicates distance to the tip
- Recommend for difficult intubation and replacing an endotracheal tube.
- Available with solid type and hollow type. Hollow type can facilitate delivery of oxygen in emergency.
- Packing:1pc/blister bag,100pcs/ctn

	Endotrachealtubus Introducer
	Endotrachéale Tube Introducteur
	Estilete de intubación
	気管内チューブイントロデューサ



OFFICE:

Jumeirah Lake Towers, Cluster G, JBC 1
Office 3801, Dubai, United Arab Emirates

info@gsticmedical.com
www.Qcaremd.com